



# MedicalDirector<sup>®</sup>

## MedicalDirector Clinical/PracSoft 3.16a Release Notes

The following information outlines the enhancements, modifications and fixes contained in the MedicalDirector Clinical/PracSoft 3.16a update, released early 2016. This information is also available for download as a PDF file from the MedicalDirector Knowledge Base at [www.medicaldirector.com/kb](http://www.medicaldirector.com/kb). The 3.16a release is an *update*, and can be installed over 3.16.





## Content for MedicalDirector Clinical

---

### Immunisations

---

- In line with amendments to the 10<sup>th</sup> Immunisation Handbook, patients with 'Chronic Liver Disease' are now identified as being at risk of influenza. The Influenza 'At Risk' Search module has been updated to include 'Chronic Liver Disease' as one of the search criteria. From the MedicalDirector Clinical Front Screen, select **Search > Influenza 'At Risk'**.

### Letter Writer

---

- Letter Writer now correctly displays Investigations Result content where the result itself contains embedded HTML or RTF information.

### Patient Details

---

- The **Alcohol** tab of the patient details window now incorporates the Audit-C assessment, for helping to identify persons who are hazardous drinkers or have active alcohol use disorders (including alcohol abuse or dependence).

### Tool Box

---


- **Modification to the Absolute CVD Risk Calculator:** For patients considered 'high risk', you can simply click the **Save** button to automatically document the risk value as >15% in the Tool Box. This will also add a note to the patient's Progress Notes.



## Content for MedicalDirector PracSoft

### Appointment Book Tools

- The tool 'Recalculate All Descriptions' has been renamed "Recalculate All Future Descriptions" to better reflect its functionality.



### Financial

- Fixed an issue when the Data Integrity Check re-calculated fees on items with GST, and generated outstanding amounts in the patient's account.
- Fixed a discrepancy where the total on the invoice did not match the total when the GST amount was added to the non-GST amount.
- Fixed an issue where a receipt for \$0.00 was printed, and remained indefinitely in Processing > Banking > Calculate list.
- Fixed an issue where, when editing an in-hospital visit that had been placed on hold, any custom Service Item Fees reverted to the default fee value.

### General

- Fixed the issue where a site experiences ongoing frequent intermittent, connectivity related issues, predominantly in the Appointment book.

MedicalDirector Head Office  
 Level 5, 72 Christie Street,  
 St Leonards NSW, 2065  
 Tel: +61 2 9906 6633  
 Fax: +61 2 9906 8910  
[www.medicaldirector.com](http://www.medicaldirector.com)  
 A.B.N. 76 068 458 515  
 Copyright© 2016