



MedicalDirector[®]
BLUECHIP

Restoring a MedicalDirector Blue Chip Database

The following outlines the process for restoring (or 'attaching') a customer's MedicalDirector Blue Chip Database, giving them the ability to log into either the MedicalDirector Blue Chip Demonstration database or their Live database.

Before you Begin

Please contact your System Administrator if you need assistance with the following:

- Ensure you have logged into Windows with local Administrative privileges.
- Ensure you have installed MedicalDirector Blue Chip on the Server computer.
- If you are using Full SQL 2005 or 2008, you will need to download the required SQL Studio Management Tools from the Microsoft web site.

Install Microsoft SQL Studio Management

You will need to be on the Server computer. The 2005 SQL Studio Management installation file is installed with a default MedicalDirector Blue Chip installation, and can be found in the following location:

- C:\Program Files\Health Communication Network\Blue Chip\Tools\
- On the MedicalDirector Blue Chip CD, browse to \Microsoft SQL Server 2005 Express\SQLServer2005_SSMSEE.msi

There are two versions (32-bit and 64-bit) of 2005 SQL Studio Management Studio Express. These are for 2005 SQL installations only - if you are using SQL 2008 please refer to **Before You Begin**, above.

- 32 bit - SQLServer2005_SSMSEE.msi
- 64 bit - SQL Server2005)SSMSEE_x64msi

Restoring the MedicalDirector Blue Chip Database

1. Browse to the location of the your MedicalDirector Blue Chip Database
2. Copy (do not 'cut') the MedicalDirector Blue Chip Database
3. Paste the MedicalDirector Blue Chip Database either of the following locations:
 - For 32-bit - C:\Program Files\Microsoft SQL Server\MSSQL.1\MSSQL\Data\
 - For 64-bit - C:\Program Files (x86)\Microsoft SQL Server\MSSQL.1\MSSQL\Data\
4. Start SQL Server Management Studio.
5. From the Windows Start menu, select **All Programs**. Within the list of programs, locate **Microsoft SQL Server** and then locate and click **SQL Server Management Studio**.
6. Once SQL Server Management Studio has opened, you will be presented with a Connect to Server window. Check the following is correct then click **Connect**.
 - Server Type: Database Engine.
 - Server name: Server Name/HCNSQL07
 - Authentication.
7. You will be presented with the **Object Explorer** screen.
8. Click on the + to expand the databases.
9. Right-click **Databases**, and select **Attach**.
10. Click **Add** on the **Attach Databases** window
11. The **Locate Database File** window will be presented. Select the MedicalDirector Blue Chip database.
12. Click Ok.
13. You will be presented with the **Attach Database** window. Click OK.

The MedicalDirector Blue Chip Database is now restored.

Configure MedicalDirector Blue Chip

To switch between the Demonstration database and the customer's Live database, follow the instructions below:

Creating the Customer Data Folder

1. Browse to the following location
 - XP - C:\Documents & Settings\All Users\Documents\Health Communication Network\Blue Chip\BlueData
 - Win 7/2008/2011 - C:\Users\Public\Health Communication Network\Blue Chip\BlueData
2. Locate and copy the entire BCDemo folder.
3. Paste the BCDemo folder into the same location.
 - a. Rename the copied BCDemo database to the customer's database name.

Configure the Database.ini file

1. Locate and open the **database.inf** file
2. Highlight and copy all text between and including the lines **[BCDemo]** and **DBSQLDatabaseImport=none**
3. Create space at the end of the file by pressing enter
4. Right click and Paste
5. Change all references of BCDemo in the copied data to the Customers Database name (as shown in bold below)
6. Save the changes

```
[BCDemo]
DBSQLUsername=sa
DBSQLPassword=9E0D2134DCEB08B4FE13DBF934AFB9DBC85D11C85D9284C9
DBSQLServer=Server Name\HCNSQL07
DBSQLDatabasePractice=BCDemo07
DBSQLDatabaseImport=none
```

```
[Customers Database Name]
DBSQLUsername=sa
DBSQLPassword=9E0D2134DCEB08B4FE13DBF934AFB9DBC85D11C85D9284C9
DBSQLServer=Server Name\HCNSQL07
DBSQLDatabasePractice=Customers Database Name
DBSQLDatabaseImport=none
```

```
[CurrentVersion]
BlueChip=2.8.2.0
```

Now when MedicalDirector Blue Chip is started, both the demonstration database and the customer database will be available for use.

MedicalDirector Head Office
Level 5, 477 Pitt Street,
HAYMARKET NSW, 2000
Tel: +61 2 9906 6633
Fax: +61 2 9906 8910
www.medicaldirector.com
A.B.N. 76 068 458 515
Copyright© 2017